

B-LIGHT 502

Upmarket chair that stands out for its innovative design. Equipped with 4 motors that control height and backrest, seat and legrest inclination. Reclining and removable armrests and headrest. Highly-stable steel structure with white LED light built in the base. White PU upholstery of high quality.

Camilla de alta gama que destaca por su diseño innovador.

Dispone de 4 motores que controlan la altura, la inclinación del respaldo, del asiento y del reposapiés. Reposabrazos y reposacabezas abatibles y extraíbles. Estructura de acero muy estable con luz LED blanca integrada en la base. Tapizado PU blanco de alta calidad.

Fauteuil haut de gamme qui se détache par son design novateur. Il est équipé de 4 moteurs qui contrôlent la hauteur et l'inclinaison du dossier, du siège et du repose-pied. Repo- se-bras et repose-tête rabattables et amovibles. Structure en acier très stable avec lumière LED blanche intégrée dans la base. Revêtement en PU blanc de haute qualité.

Hochwertige Liege, die durch ihr innovatives Design hervorsteicht. Diese Liege ist mit 4 Motoren ausgestattet, um die Höhe, die Rückenlehne, die Sitzneigung und Fußteil zu regulieren. Die Armlehne, Kopfstütze und Fußteil sind klappbar und abnehmbar. Sehr stabile Stahlkonstruktion mit LED weißer Beleuchtung am Gestell der Liege. Hochwertige PU Polsterung.

Letino di alta gamma, che si distingue per il suo design innovativo. Dispone di 4 motori che controllano l'altezza, l'inclinazione dello schienale, la seduta e il poggiatesta. Braccioli, poggiatesta e poggiatesta pieghevoli e removibili. Struttura in acciaio molto stabile con luce LED bianca integrata nella base. Rivestimento in PU bianco di alta qualità.

Electric Beauty Bed Ref. BLIGHT502.4.A26

Size without armrest:
193 x 56 x 60/105 cm
Size with armrest:
193 x 82,5 x 60/105 cm
Weight: 94 Kg
Upholstery: PU
Thickness mattress: 10 cm
Density mattress: 30 Kg/m³
Motors: 4

- 1 Breather hole 13,5 x 9 cm
- 2 Height (up & down) 60/92 cm
- 3 Adjustable Armrest
- 4 Backrest angle 0° ~ 60°
- 5 Legrest angle 0° ~ 90°
- 6 Trendelenburg 0° ~ 11°
- 7 White LEDs
- 8 Hand Controller Ref. HC004.WKE015

OPTIONS

Cover Weelko Beauty
Foot controller Ref. FC004.WKE015
Wheels No
Paper roll holder Ref. 50080.1

COLOURS



Packaging size cm 139 x 63 x
N.W. kg 59.94 kg
G.W. kg 105.5 kg

Loading capacity 175 kg



B-LIGHT 502

