

KUNE

Functional 3-section pedicure chair of great stability equipped with 2 motors that control height and backrest inclination. Removable and reclining armrests. Legrests swivel up and are independently adjustable in height by gas lift. Adjustable headrest cushion with breather hole. White PU upholstery of high quality that eases cleaning.

Sillón de pedicura muy estable equipado con 2 motores que controlan la altura y la inclinación del respaldo. Reposabrazos abatibles y extraíbles. Reposapiés con apertura lateral y ajustables en altura independientemente mediante pistón de gas. Con agujero facial y cojín reposacabezas regulable. Tapizado blanco en PU de alta calidad y fácil de limpiar.

Fauteuil de pédicure à 3 plans de design fonctionnel avec une structure très stable. Il est équipé de 2 moteurs qui contrôlent la hauteur et l'inclinaison du dossier. Accoudoirs amovibles et inclinables. Hauteur des jambières réglable de façon indépendante par vérin à gaz et écartables. Coussin repose-tête ajustable avec trou visage. Revêtement en PU blanc de haute qualité et facile à nettoyer.

Stabile Fußpflegeliege mit 2 Motoren ausgestattet, um die Höhe und die Neigung der Rückenlehne zu regulieren. Abklappbare und abnehmbare Armlehne. Die Fußstützen können mittels Gasfeder und unabhängig voneinander in der Höhe verstellt werden. Sie lassen sich nach Aussen drehen. Weiße PU-Polsterung, leicht abwaschbar.

Poltrona per pedicura molto stabile, dotata di 2 motori che controllano l'altezza e l'inclinazione dello schienale. Braccioli pieghevoli e removibili. Poggiatesta con apertura laterale e regolabili in altezza in maniera indipendente tramite pistone a gas. Con foro facciale e cuscino poggiatesta regolabile. Rivestimento in PU bianco di alta qualità e facile da pulire.

Electric Podiatry Chair Ref. WKP007.2.A26

Size with armrest: 185x79
Size without armrest: 185x58
Height: 58/81cm
Weight: 51 Kg
Upholstery: PU
Thickness mattress: 8 cm
Motors: 2

MULTIPLE
LEG
POSITIONS



8

- 1 Breather hole + cushion
Breather hole 14 x 12 cm
- 2 Height (up & down) 58/81 cm
- 3 Armrest rotation 0-180°
- 4 Backrest angle 0° ~ 70°
- 5 Legrest angle 0° ~ 90°
- 6 Leg extension ±16 cm
- 7 Hand controller Ref. HC002.WKP007

OPTIONS

Cover Weelko Podo
Foot controller Ref. FC002.WK
Wheels No
8 Paper roll holder Ref. 50080.1

COLOURS



Packaging size cm 152 x 70 x 64
N.W. kg 51 kg
G.W. kg 62kg

Loading capacity 150 kg



KUNE

